

Catalog of medical consumables



Catalogs of instruments



You can request these catalogs from our sales representatives or order it on e-mail address medilab@medilab.cz

- Catalog of basic surgical instruments
- Catalog of orthopedic instruments
- Catalog of instruments for gynecology

Index

Syringes

Hypodermal needles

Infusion sets

Infusion set

Infusion set with port

Extension tubes

Extension tubes

2-way extension tubes

Transfusion set

Peripheral intravenous catheters

Peripheral intravenous catheters — butterfly type

Peripheral intravenous catheters — pen type

Safety connectors and stopcocks

Heparin caps

3-way stopcocks

Flow-in — safety injection site adapter

Urinal bags

Urinal bags

Urine bag holders

Urinal bag — leg type

Fastening stripes for urinal bags — leg type

Foley catheters

2-way Foley catheters Nelaton for women

2-way Foley catheters Nelaton for children

Closure for Foley catheter

2-way Foley catheters Nelaton for men and women

2-way Foley catheters Tiemann

Oxygen masks

Oxygen nasal cannulas and airways

Bedding and sheets

Disposable blanket

Disposable pillow

Disposable nonwoven sheets

Bed sheets and bed covers

Disposable bed sheets — nonwoven textile

Gloves and medical brush

Latex examination gloves

Surgical cleaning brush

Surgical face masks

Surgical caps

Surgical clothing

Surgical shirts

Surgical pants

Cytology brushes

Laboratory equipment

Microscopic slides

Histological material containers

Biological material containers

Pipettes

Microtubes

Tubes

Polyethylene lid

Tubes with conical bottom

Microscope slide boxes

Cardboard slide mailers

Archiving box for histological cassettes

Container tray

HistoMED

Pre-filled tubes HistoMED

Pre-filled containers HistoMED





Surgical suture material

Syringes




3-part syringes

Disposable 3-part syringes are packed in practical sterile peel-bag. They are equipped with double gasket, which prevents air penetration and ensures fluent move of plunger in the syringe cylinder. Syringes are intended for dosage of medicals, collection of body fluids and other medical acts.



Product description	REF	pcs/box	pcs/carton
 Insulin syringe U100 – 0,5 ml	IS0529G	100	3 600
 Insulin syringe U100 – 1 ml	IS1029G	100	3 600
 1 ml, sterile, centric cone, Luer-slip, non-residual	KD1301B	100	3 600
 2 ml, sterile, centric cone, Luer-slip	KD1302	100	2 700

Syringes



3-part syringes

Disposable 3-part syringes are packed in practical sterile peel-bag. They are equipped with double gasket, which prevents air penetration and ensures fluent move of plunger in the syringe cylinder. Syringes are intended for dosage of medicals, collection of body fluids and other medical acts.



Product description	REF	pcs/box	pcs/carton
 50 ml, sterile, centric cone, Luer-lock	KD1352	50	300
 60 ml, sterile, centric cone, with cone	MA1450P	25	200

Hypodermal needles



Hypodermal needles

Disposable hypodermal needles are color-coded according to international codes. Needles are packed in sterile peel-bag.


Product description			REF	pcs/box	pcs/carton
 25G x 1"	0,50 mm x 25 mm	orange	HN25G	100	10 000
 23G x 1"	0,60 mm x 25 mm	blue	HN23G	100	10 000

Infusion sets



Infusion sets (length 150 cm)

Infusion sets are intended for administration of infusion solutions from bottles or plastic pouches into the patient's vein by gravity. Infusion sets are packed in single sterile package and are intended for single use. Without phthalates.


Product description	REF	pcs/box	pcs/carton
 Infusion set for gravity infusion	MA2208P	25	400

Infusion sets




Infusion set with port (length 150 cm)

Infusion sets are intended for administration of infusion solutions from bottles or plastic pouches into the patient's vein by gravity. Infusion sets are packed in single sterile package and are intended for single use. Without phthalates.





Product description	REF	pcs/box	pcs/carton
 Infusion set for gravity infusion with port	MA1108PORT	25	400

Extension tubes




Extension tubes

Extension tubes can be used in combination with infusion sets, transfusion sets and peripheral intravenous catheters. Extension tube SH10FLW is with safety needle free connector.


Product description	REF	pcs/box	pcs/carton
 Extension tubes 1,8 × 450	SH45	100	800
 Extension tubes 1,8 × 100	SH10FLW	50	600
 Extension tubes 1,8 × 1800	SH180	50	400
 2-way extension tubes	SHY	100	800

Transfusion set



Transfusion set

Transfusion set is intended for transfer of stored blood, primary blood components, blood substitutes etc. from bottles or pouches into patient's vein. Without phthalates.

Product description	REF	pcs/box	pcs/carton
 Transfusion set	KD301N	25	250

Peripheral intravenous catheters



Peripheral intravenous catheters — butterfly type

Intravenous catheters — with flexible wings for easy fixation, luer-lock connector, integrated three radiocontrast stripes, transparent chamber for blood return control, atraumatic tip shape. PVC free.




Product description	REF	pcs/box	pcs/carton
● 18G — green — butterfly type, dimensions 1,2 × 45 mm	PZK18GB	50	1 000
● 20G — pink — butterfly type, dimensions 1,0 × 33 mm	PZK20GB	50	1 000
● 22G — blue — butterfly type, dimensions 0,8 × 25 mm	PZK22GB	50	1 000
● 24G — yellow — butterfly type, dimensions 0,7 × 19 mm	PZK24GB	50	1 000

Peripheral intravenous catheters



Peripheral intravenous catheters — „pen type“

Intravenous catheters — without wings, without port — „pen type“ — luer-lock connector, integrated three radiocontrast stripes, transparent chamber for blood return control, atraumatic tip shape. PVC free.

Product description	REF	pcs/box	pcs/carton
 20G — pink — pen type, dimensions 1,0 x 33 mm	PZK20PL	50	2 000
 22G — blue — pen type, dimensions 0,8 x 25 mm	PZK22PL	50	2 000
 24G — yellow — pen type, dimensions 0,7 x 19 mm	PZK24PL	50	2 000

Safety connectors and stopcocks



Heparin caps




Heparin cap can be connected to peripheral venous catheter or connection tube for continuous distribution of medication.

3-way stopcocks

Multi-way ventile for infusion sets and monitoring, lipidoresistant.

Flow-in — safety injection site adapter

Automatic ventile for closure of peripheral intravenous catheters and connection tubes.



Product description	REF	pcs/box	pcs/carton
 Heparin caps	HCY	100	5 000
 3-way stopcocks	STF11	50	1 000
 Flow-in	FLW	100	2 000

Urinal bags



Urinal bags

Urinal bags are suitable for collection of secrets, capacity 2000 ml, sterile with stop-return valve and length of tube 90 or 120 cm. Bags can be hung on patient's bed. Urinal bag without outlet must be emptied after reaching the bag capacity and discarded on determined place. Urinal bag with outlet is intended for repeatable usage for one patient, after usage must be emptied and discarded on determined place. All bag types have suspension holes for hanging the bag with the bag holder on the patient's bed.

Product description	REF	pcs/box	pcs/carton
 Urinal bags with bottom cross valve, 2000 ml, length of tube 90 cm, sterile single packing	UBT2000	10	250
 Urinal bags with bottom cross valve, 2000 ml, length of tube 120 cm, sterile single packing	UBT2000L	10	250
 Urinal bags without outlet, 2000 ml, length of tube 90 cm, sterile single packing	UBT2000B	10	250
 Urine bag holder - hanger	UPHR	50	500
 Urine bag holder	UBH	100	2 000

Urinal bags





Urinal bag – leg type

Urinal bags are suitable for collection of secrets, capacity 750 ml, three chambers, with stop-return valve and length of tube 45 cm, fulfills the needs of every-day usage. Bags are fastened with re-usable set of two velcro stripes, which is not included in the packing. Rear part of the bag is textile coated for sweat absorption and skin comfort. Urine is distributed in three chambers, the bag remains flat for maximum discretion. Distribution of content in chambers helps to safe fastening during moving.

Fastening stripes for urinal bags – leg type

Set of two velcro stripes for fastening the urinal bag to patient's leg. Stripes are coated with silicone stripe, which prevents slipping off. Fastening stripes are made of flexible latex-free material.

Product description	REF	pcs/box	pcs/carton
 Urinal bag – leg type, capacity 750 ml, three chambers, length of tube 45 cm, single packing	ULB500	10	200
 Set of fastening velcro stripes for urinal bags – leg type, length 50 cm, width 2,5 cm, white	LBS	10 pairs	500 pairs

Foley catheters




2-way Foley catheters Nelaton for women (length 20 cm)

Permanent balloon catheter is intended for catheterisation of urine bladder. Sterile, 100% silicone coated.

Product description	REF	pcs/box	pcs/carton
○ Size CH12, balloon capacity 5 – 10 ml, for women	2FC12FRW	10	500
● Size CH14, balloon capacity 5 – 10 ml, for women	2FC14FRW	10	500
● Size CH16, balloon capacity 5 – 15 ml, for women	2FC16FRW	10	500
● Size CH18, balloon capacity 5 – 15 ml, for women	2FC18FRW	10	500
● Size CH20, balloon capacity 5 – 15 ml, for women	2FC20FRW	10	500
● Size CH22, balloon capacity 5 – 15 ml, for women	2FC22FRW	10	500

Foley catheters



2-way Foley catheters Nelaton for children (length 20 cm)


Permanent balloon catheter is intended for catheterisation of urine bladder. Sterile, 100% silicone coated.

Closure for Foley catheter „spigot“

Intended for closing of foley catheter (Nelaton, Tiemann).

Product description	REF	pcs/box	pcs/carton
● Size CH8, balloon capacity 3 ml, for children	2FC8FRM3	10	500
● Size CH10, balloon capacity 3 ml, for children	2FC10FRM3	10	500
● Closure for Foley catheter - „spigot“ sterile	UFK	100	1 000

Foley catheters




2-way Foley catheters Nelaton (long, length 40 cm, for men and women)

Permanent balloon catheter for catheterization of urine bladder. Sterile, 100% silicone coated.

Product description	REF	pcs/box	pcs/carton
○ Size CH12, balloon capacity 5 – 15 ml, for men and women	2FC12FRM5	10	500
● Size CH14, balloon capacity 5 – 15 ml, for men and women	2FC14FRM5	10	500
● Size CH16, balloon capacity 5 – 15 ml, for men and women	2FC16FRM5	10	500
● Size CH18, balloon capacity 5 – 15 ml, for men and women	2FC18FRM5	10	500
● Size CH20, balloon capacity 5 – 15 ml, for men and women	2FC20FRM5	10	500
● Size CH22, balloon capacity 5 – 15 ml, for men and women	2FC22FRM5	10	500
● Size CH24, balloon capacity 15 – 30 ml, for men and women	2FC24FRM5	10	500
● Size CH26, balloon capacity 15 – 30 ml, for men and women	2FC26FRM5	10	500

Foley catheters



2-way Foley catheters Tiemann (long, length 40 cm)

Permanent balloon catheter for catheterization of urine bladder. Sterile, 100% silicone coated.

Product description	REF	pcs/box	pcs/carton
○ Size CH12, balloon capacity 5 – 15 ml	2FC12TM	10	500
● Size CH14, balloon capacity 5 – 15 ml	2FC14TM	10	500
● Size CH16, balloon capacity 5 – 15 ml	2FC16TM	10	500
● Size CH18, balloon capacity 5 – 15 ml	2FC18TM	10	500
● Size CH20, balloon capacity 5 – 15 ml	2FC20TM	10	500
● Size CH22, balloon capacity 5 – 15 ml	2FC22TM	10	500
● Size CH24, balloon capacity 5 – 15 ml	2FC24TM	10	500

Oxygen masks



Oxygen masks — standard / nonrebreathable / with nebulizer

Oxygen masks are intended for administration of oxygen therapy. Sterile, with tube length 2 m.

Product description	REF	pcs/package	pcs/carton
 Oxygen mask — standard — size L, for adolescents	KMSL	1	100
 Oxygen mask — standard — size XL, for adults	KMSXL	1	100
 Oxygen masks — nonrebreathable, size L, for adolescents	KMVL	1	50
 Oxygen masks — nonrebreathable, size XL, for adults	KMVXL	1	50
 Oxygen mask — nebulizer — size L, for adolescents	KMRL	1	50
 Oxygen mask — nebulizer — size XL, for adults	KMRXL	1	50

Oxygen nasal cannulas and airways





Oxygen nasal cannulas

Oxygen nasal cannulas for adults. Tube is equipped with cross-section, which prevents from its complete cranking. Length 2 m.

Airways

Airways are anatomically shaped, compact, disposable and sterile. Different colors for size distinction.

Product description	REF	pcs/package	pcs/carton
 Oxygen nasal cannulas for adults — size L	KB	20	200
 Airways — size 3, length 90 mm, yellow	GA3	1	50
 Airways — size 4, length 100 mm, red	GA4	1	50

Bedding and sheets



Disposable blanket

Disposable blanket with polyester filling is made of non-woven fabric. The blanket is waterproof.

Disposable pillow


Disposable pillow with polyester filling is made of non-woven fabric. It is suitable for ARO, ICU, infection department or rescue services.

Disposable nonwoven sheets

The sheets are made of non-woven fabric. They are intended to use as a one-time hygienic protection for hospital bedding.

Product description	REF	pcs/box	pcs/carton
 Disposable blanket 190 × 110 cm, 400 g, blue	DJ400B	1	48
 Disposable nonwoven blanket cover, 195 × 115, 25 g	JP25195	10	100
 Disposable pillow 40 × 40 cm, white	PJ40x40-W	1	36
 Disposable pillow-case, 42 × 42 cm, 25 g	JP2542	10	200

Bed sheets and bed covers



Disposable bed sheets - nonwoven textile

Bedsheets can be used as disposable hygienic protection on hospital and examination beds. Bedsheet REF PRO3015200 is additionally equipped with polyethylene waterproof layer.

Product description	REF	pcs/package	pcs/carton
Disposable bed sheet, nonwoven textile, 210 × 90 cm, 60 g/m ²	PRO60210	10	100
Disposable bed sheet, nonwoven textile, 200 × 80 cm, 40 g/m ²	PRO40200	10	100
Disposable bed sheet, nonwoven textile, 220 × 80 cm, 25 g/m ²	PRO25220	10	100
Disposable bed sheet, nonwoven textile, 200 × 80 cm, PP 30 g/m ² + PE 15 g/m ²	PRO3015200	10	100

Gloves and medical brush



Latex examination gloves

Examination gloves are left-right and they are for single use. They distinguish in high durability as well as in required softness. They are used for examination, nursing, diagnostic, therapeutic and other procedures in health sector.

Surgical cleaning brush

Soft bristles, which reliably clean joints, nails and skin, thus being really gentle. Brush has integrated sponge, which reliably cleans hands and forearms without danger of irritation or damage (microtrauma) of skin. Nail cleaner included.

Product description	REF	pcs/box	pcs/carton
 Latex examination gloves — size S, powdered	LEGSPP	100	1 000
 Latex examination gloves — size M, powdered	LEGMPP	100	1 000
 Latex examination gloves — size L, powdered	LEGLPP	100	1 000
 Latex examination gloves — size S, powder-free	LEGSPF	100	1 000
 Latex examination gloves — size M, powder-free	LEGMPF	100	1 000
 Latex examination gloves — size L, powder-free	LEGLPF	100	1 000
 Surgical cleaning brush with sponge and nail cleaner, sterile	MB	30	300

Surgical face masks





Surgical face masks

Bacterial efficiency of surgical face masks is $\geq 98\%$. Material is non-irritating and skin-comfortable. Face mask with ties is fastened with binding ties and is equipped with adjustable strip in the nose and mouth area.

Face mask with elastic strips is fastened with comfortable elastic strips, this face mask is also equipped with adjusting strip in the nose and mouth area.

All face masks comply with standard EN 14683. Protection type II.

Product description	REF	pcs/box	pcs/carton
 3-ply face mask with ties, blue	CHOM1	50	2 000
 3-ply face mask with elastic strips, blue	CHOM2	50	2 000

Surgical face masks



Surgical face masks – with flowers / with hearts / black with silver flowers

Bacterial efficiency of surgical face masks is $\geq 98\%$. Material is non-irritating and skin-comfortable.

Face mask with ties is fastened with binding ties and is equipped with adjustable strip in the nose and mouth area.

Face mask with ear-loops is fastened with comfortable elastic strips, this face mask is also equipped with adjusting strip in the nose and mouth area.

All face masks comply with standard EN 14683. Protection type II.

Product description	REF	pcs/box	pcs/carton
 3-ply face mask with ties, blue flowers	CHOM5	50	2 000
 3-ply face mask with ear-loops, blue flowers	CHOM6	50	2 000
 3-ply face mask with ties, pink flowers	CHOM7	50	2 000
 3-ply face mask with ear-loops, pink flowers	CHOM8	50	2 000
 3-ply face mask with ties, pink hearts	CHOM9	50	2 000
 3-ply face mask with ear-loops, pink hearts	CHOM10	50	2 000
 3-ply face mask with ties, black with silver flowers	CHOM11	50	2 000
 3-ply face mask with ear-loops, black with silver flowers	CHOM12	50	2 000

Surgical caps



Surgical cap – with ties




Disposable blue surgical cap with ties is made of soft breathable nonwoven textile. Surgical cap is suitable for using during surgical operation as well as universal head cover.

Surgical cap – with ties and absorption strip

Disposable blue surgical cap with ties is made of soft breathable nonwoven textile. Surgical cap is equipped with absorption stripe. It is suitable for using during surgical operation as well as universal head cover.

Surgical cap – with elastic strip

Disposable blue surgical cap is made of soft breathable nonwoven textile. Surgical cap has elastic strip in the rear part and is suitable for using during surgical operation as well as universal head cover.

Product description	REF	pcs/box	pcs/carton
 Surgical cap with ties, blue	SCTOB	100	1 000
 Surgical cap with ties and absorption strip, blue	SCWAS	100	1 000
 Surgical cap with elastic strip, blue	SCPB	100	1 000

Surgical caps



Surgical caps — with flowers / with hearts / black with silver flowers

Disposable surgical caps are available with ties or elastic strip. They are made of soft breathable nonwoven textile. Surgical caps are suitable for using during surgical operation as well as universal head cover.

Product description	REF	pcs/box	pcs/carton
 Surgical cap with ties, blue flowers	SCTOPWB	100	1 000
 Surgical cap with ties, pink flowers	SCTOPWP	100	1 000
 Surgical cap with ties, pink hearts	SCTOPH	100	1 000
 Surgical cap with ties, black with silver flowers	SCTOPBS	100	1 000
 Surgical cap with elastic strip, blue flowers	SCPWB	100	1 000
 Surgical cap with elastic strip, pink flowers	SCPWP	100	1 000
 Surgical cap with elastic strip, pink hearts	SCPH	100	1 000
 Surgical cap with elastic strip, black with silver flowers	SCPBS	100	1 000

Surgical clothing

















Surgical shirt — flowered pocket / pocket with hearts

Disposable surgical shirts are made of soft, breathable non-woven fabric. Shirts are suitable to be worn during surgery. They are available in sizes S, M, L, XL and variants with a flowered pocket or a pocket with hearts. The shirts have short sleeves and a V-neck.


Surgical pants

Disposable surgical pants are straight-cut with a tie at the waist, are made of soft, breathable non-woven fabric. Pants are suitable to be worn during surgery and available in sizes S, M, L and XL.





Product description	REF	pcs/box	pcs/carton
 Surgical shirt — pink flowers S	OSSPF		50
 Surgical shirt — pink flowers M	OSMPF		50
 Surgical shirt — pink flowers L	OSLPF		50
 Surgical shirt — pink flowers XL	OSXLPF		50
 Surgical shirt — blue flowers S	OSSBF		50
 Surgical shirt — blue flowers M	OSLBF		50
 Surgical shirt — blue flowers L	OSMBF		50
 Surgical shirt — blue flowers XL	OSXLBF		50
 Surgical shirt — pink hearts S	OSSH		50
 Surgical shirt — pink hearts M	OSMH		50
 Surgical shirt — pink hearts L	OSLH		50
 Surgical shirt — pink hearts XL	OSXLH		50
 Surgical pants — blue S	OTS		50
 Surgical pants — blue M	OTM		50
 Surgical pants — blue L	OTL		50
 Surgical pants — blue XL	OTXL		50

Cytology brushes



Cytology brushes

Soft sample collection brush on plastic shaft, length 20 cm.

Product description	REF	pcs/box	pcs/carton
 Cytology brush — sterile	CBS	100	5 000
 Cytology brush — nonsterile	CB	100	5 000

Laboratory equipment



Microscope slides






Microscope slides with both sides rounded, clear or one side frosted. 45° tapered corners.

Microscope slides — printable, positive charged

Microscope slide — resistant to commonly used laboratory solvents. 45° tapered corners.

Positive charged microscope slides — these slides have permanent positive charge. The static electricity pulls parts of frozen tissues and cytology preparations on slide surface. Parts of tissues and cytology preparations adhere better to positive charged slides, without need of special glue or protein coating. Adhesive variant of microscope slide covered by silane.

*) Covering glass - vacuum packed 500 pcs.

Product description	REF	pcs/box	pcs/carton
 Microscope slides — clear, corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	ME7101	50	2 500
 Microscope slides — one side frosted, corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	ME7105-1	50	2 500
 Microscope slides, printable — white, corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	ME7109-W	50	2 500
 Microscope slides, printable — white, corners 90°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	ME7109-W90	50	2 500
 Microscope slides — positive charged — green, corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	MEPN-G	50	2 500
 Microscope slides — positive charged — pink, corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	MEPN-P	50	2 500
 Microscope slides — positive charged — yellow, corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	MEPN-Y	50	2 500
 Adhesive microscope slides, covered by silane, positive charged — blue, corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	MEAS-B	50	2 500
 Covering glass, dimensions 24 × 50 mm, thickness 0,13 — 0,17 mm	MCG7201	100 *)	20 000

Laboratory equipment












Histological material containers


Containers for histological material are made of nontoxic polypropylene. Cap is made of HDPE. Suitable for collection, transport and storage of liquid samples. Prevents leakage and vaporization of formalin. Upper and bottom part ribbed for easy opening and closing with gloves. Scales on three sides.

Biological material containers

Containers for biological material are made of virgin polystyrene and polypropylene. All materials are non-toxic. Size 7 ml is suitable for collection of small samples. Size 30 ml has unique cap with protection against leakage.

Product description	REF	pcs/carton
 Histological material container, capacity 20 ml	28100	1 000
 Histological material container, capacity 40 ml	28200	500
 Histological material container, capacity 60 ml	28300	500
 Histological material container, capacity 90 ml	28400	300
 Histological material container, capacity 120 ml	28500	300
 Histological material container, capacity 160 ml	28900	300
 Histological material container, capacity 250 ml	28600	154
 Histological material container, capacity 500 ml	28700	132
 Histological material container, capacity 1000 ml	28800	64
 Biological material container, capacity 7 ml	22500	700
 Biological material container, capacity 30 ml	22100	400

Laboratory equipment



Pipettes

Disposable pipette is made of special nontoxic low density polyethylene. It is suitable for safe transfer and distribution of liquids in all laboratory types. Eliminated risk of infection transfer.

Microtubes








Microtubes are made of robust virgin polypropylene, they are transparent for better visibility. Tubes with safety lock have perfect stability, tightness, protection against vaporization and unwanted opening. Flat lid with tapered edges enables easy opening and closing. Frosted part for writing and marking. Moulded transparent scale. Lid with writable surface. Sterilizable in autoclave.

Tubes

Tube is capable of centrifugation speed over 3000 x g and temperature to 121°C.

Polyethylene lid

The polyethylene lids are stoppers with hollow upper part with adjustment for gripping. They fit in test tubes and also in round cuvettes and centrifugation tubes. 3 inner ribs for safe fastening to tubes.

Product description	REF	pcs/package	pcs/carton
 Disposable pipette, 4 ml	13301	500	3 000
 Microtube – 1,5 ml, without lid	36210	1 000	11 000
 Microtube with scale, with lid	35200N	500	10 000
 Tube 16 x 100 mm, volume 14 ml, without lid	81621	500	2 500
 Polyethylene lid for tube, red and blue color	89216R	1 000	16 000
 Tube with conical bottom, free standing, volume 50ml, with lid	HP1001	25	500
 Tube with conical bottom, volume 50ml, with lid	HP1003	25	500

Laboratory equipment



Microscope slide boxes

Plastic boxes are used to store and transfer microscope slides and prevent them from breaking or sticking.

Cardboard slide mailers


Slide mailers are made of cardboard and designed for holding microscope slides.











Archiving box for histological cassettes

Cardboard box is designed for storing histological cassettes.

Container tray

The polystyrene tray is used to store containers with a diameter of 45 and 32 mm.



Product description	REF	pcs/package	pcs/carton
 Microscope slide box for 5 pcs of slides, transparent	HP4007	4	100
 Microscope slide box for 12 pcs of slides (orange colour)	HP4008	2	
 Microscope slide box for 50 pcs of slides	BFCS	1	100
 Microscope slide box for 100 pcs of slides	BFCS-L	10	50
 Slide mailer 295 x 205 mm	DEKLS	10	40
 Slide mailer 230 x 460 mm	DEKLM	10	30
 Slide mailer 230 x 485 mm	DEKLL	10	30
 Slide mailer — pad without lid	DEKL-BV	10	40
 Archiving box for histological cassettes 300 x 420 x 50 mm	AKE300x420x50		50
 Polystyrene tray for containers with diameter 45 and 32 mm	XPS255x205x50	1	

HistoMED



Pre-filled tubes and containers HistoMED

HistoMED pre-filled containers are intended for fixation and transport of samples (small excisions and biopsies) for histopathological examination. The container is always half filled with fixative solution. The container contains max. 4 % formaldehyde buffered 10% formalin, pH 7,2 - 7,4.

Product description	REF	pcs/package	pcs/carton
 Pre-filled tube HistoMED, volume 7 ml, amount of formalin 3,5 ml	HMED07	50	1 000
 Pre-filled tube HistoMED, volume 14 ml, amount of formalin 7 ml	HMED14	100	800
 Pre-filled container HistoMED, volume 20 ml, amount of formalin 10 ml	HMED20	60	120
 Pre-filled container HistoMED, volume 40 ml, amount of formalin 20 ml	HMED40	35	140
 Pre-filled container HistoMED, volume 60 ml, amount of formalin 30 ml	HMED60	35	105
 Pre-filled container HistoMED, volume 90 ml, amount of formalin 45 ml	HMED90	35	70
 Pre-filled container HistoMED, volume 120 ml, amount of formalin 60 ml	HMED120	23	46
 Pre-filled container HistoMED, volume 160 ml, amount of formalin 80 ml	HMED160	23	23
 Pre-filled container HistoMED, volume 250 ml, amount of formalin 125 ml	HMED250	9	27
 Pre-filled container HistoMED, volume 500 ml, amount of formalin 250 ml	HMED500		10
 Pre-filled container HistoMED, volume 1000 ml, amount of formalin 500 ml	HMED1000		5
 Pre-filled container HistoMED, volume 3000 ml, amount of formalin 1500 ml	HMED3000	1	8
 Pre-filled container HistoMED, volume 5000 ml, amount of formalin 2500 ml	HMED5000	1	4

Absorbable surgical suture material



	<p>Surgicryl® 910 Polyglactine 910</p> <p>Support of wound approximately 30 days</p>		
	<p>Surgicryl® PGA Polyglycolic acid</p> <p>Support of wound approximately 30 days</p>		
	<p>Surgicryl® Rapid Polyglycolic acid</p> <p>Support of wound approximately 10 days</p>		
	<p>Surgicryl® Monofast Polyglecaprone 25</p> <p>Support of wound approximately 20 days</p>		
	<p>Surgicryl® Monofilament Polydioxanone</p> <p>Support of wound approximately 60 days</p>		

Non-absorbable surgical suture material



You can request the catalog from our sales representatives or order by e-mail medilab@medilab.cz.
For more information, visit www.medilab.cz.

Certification

Certifikát

TÜVRheinland

Certifikační orgán
TÜV Rheinland LGA Products GmbH

Titulek potvrzení, že organizace
MEDILAB ČR s.r.o.
Záhumná 574
51771 Česká Meziříčí
Česká Republika

Má zaveden a uplatňuje systém řízení jakosti ve výrobě zdravotnických prostředků v následujícím rozsahu:

Výroba a distribuce sterilních a nesterilních zdravotnických prostředků na jedno použití

Byly předloženy důkazy o tom, že jsou plněny požadavky normy

EN ISO 13485:2016

Systém řízení jakosti je předmiotem každoročního dohledového auditu.

Datum účinnosti: 2019-11-22
Registrační číslo certifikátu: SX 60144023 0001
Číslo zprávy o provedeném auditu: 26300449 002
Platnost tohoto certifikátu do: 2022-04-08

Daks
Doklad
Attestace
Akreditační orgán
0-30-14189-0-02

Certifikační orgán
D. Swiatko

TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
Tel.: +49 201 800 121 Fax: +49 201 800 1000 email: service@lga.tuev-rheinland.com

Certifikát ES

Směrnice 93/42/EHS Příloha V

Zajištění jakosti výroby zdravotnických prostředků

TÜVRheinland

Registrační číslo: DD 60144023 0001
Zpráva z auditu číslo: 26300449 002

Výrobce: **MEDILAB ČR s.r.o.**
Záhumná 574
51771 Česká Meziříčí
Česká Republika

Výrobky: (zařazené výrobky jsou uvedeny v příloze)

Datum vypršení platnosti: 2024-04-08
Notifikovaný orgán tímto prohlašuje, že pro uvedené výrobky byly splněny požadavky přílohy V směrnice 93/42/EHS. Výše uvedený výrobce má zaveden a uplatňuje systém zajištění jakosti, který je objektivně prováděnatelný dohledu v souladu s přílohou V, oddílem 4 výše citované evropské směrnice. Pro provedení na trh se pro výrobky třídy II a III, které jsou předmětem tohoto certifikátu, požaduje certifikát o splněním provedení evropské úpravy šrotků podle přílohy II směrnice.

Datum účinnosti: 2019-11-22
Datum vystavení: 2019-11-22

TÜVRheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg

TÜV Rheinland LGA Products GmbH je Notifikovaný orgán podle směrnice 93/42/EHS a zkratkových pravidlech. Jeho identifikační číslo Notifikovaného orgánu je 0197.

D. Swiatko

TÜV Rheinland

LGA Products GmbH

Tillystraße 2, 90431 Nürnberg

Doc. 1/1, Rev. 0

Příloha k certifikátu

s registračním číslem: DD 60144023 0001
Zpráva z auditu č.: 26300449 002

Výrobce: **MEDILAB ČR s.r.o.**
Záhumná 574
51771 Česká Meziříčí
Česká Republika

Zařazené výrobky:

- Sterilní bodupouštní usávky Flow-in
- Sterilní 3-cestné uzavírací koberky
- Sterilní foliové katetry BELLON
- Sterilní insulínové stříkačky
- Sterilní injekční stříkačky
- Sterilní katetry jedlé
- Sterilní tenzomoduly
- Sterilní používání modulů
- Sterilní foliové katetry TIBBON

Pro následující zdravotnické prostředky potvrzuje certifikace posouzení aspekty výroby spojené se zajištěním a udržením sterility prostředků:

- Chirurgické kartáčky
- Určité sáčky

Datum: 2019-11-22

D. Swiatko

Certificate

TÜVRheinland

The Certification Body of
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization
MEDILAB ČR s.r.o.
Záhumná 574
51771 Česká Meziříčí
Česká republika

has established and applies a quality management system for medical devices for the following scope:

Production and distribution of sterile and non-sterile disposable medical devices

Proof has been furnished that the requirements specified in

EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2019-11-22
Certificate Registration No.: SX 60144023 0001
An audit was performed. Report No.: 26300449 002
This Certificate is valid until: 2022-04-08

Daks
Doklad
Attestace
Akreditační orgán
0-30-14189-0-02

Certification Body
D. Swiatko

TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
Tel.: +49 201 800 121 Fax: +49 201 800 1000 email: service@lga.tuev-rheinland.com

EC Certificate

Directive 93/42/EEC Annex V

Production Quality Assurance

Medical Devices

TÜVRheinland

Registration No.: DD 60144023 0001
Report No.: 26300449 002

Manufacturer: **MEDILAB ČR s.r.o.**
Záhumná 574
51771 Česká Meziříčí
Česká republika

Products: (see attachment for products included)

Expiry Date: 2024-04-08

The Notified Body hereby declares that the requirements of Annex V of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex V, section 4 of the aforementioned directive. For placing on the market of class IIb and class III devices covered by this certificate an EC type-examination certificate according to Annex III is required.

Effective Date: 2019-11-22
Date: 2019-11-22

TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg

TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.

D. Swiatko

TÜV Rheinland

LGA Products GmbH

Tillystraße 2, 90431 Nürnberg

Doc. 1/1, Rev. 0

Attachment to Certificate

Registration No.: DD 60144023 0001
Report No.: 26300449 002

Manufacturer: **MEDILAB ČR s.r.o.**
Záhumná 574
51771 Česká Meziříčí
Česká republika

Products included:

- Sterile injection site adapters Flow-in
- Sterile 3-way stopcocks
- Sterile Foley catheters BELLON
- Sterile insulin syringes
- Sterile syringes for single use
- Sterile infusion sets
- Sterile transfusion sets
- Sterile extension tubes
- Sterile Foley catheter TIBBON

For the following medical devices the scope covers only the aspects of manufacture concerned with securing and maintaining sterile conditions:

- Medical brushes
- Urine bags

Date: 2019-11-22

D. Swiatko

Certificates of quality management system ČSN EN ISO 13485:2016.

Certificates of product quality according to 93/42 EEC, justification for product CE marking.





MEDILAB ČR, s.r.o., Czech Republic
E-mail: medilab@medilab.cz, objednavky@medilab.cz
Phone: (+420) 494 661 067, Fax: (+420) 494 661 068

www.medilab.cz